



Recombinant Rat Prostacyclin synthase (Ptgis), partial, Yeast

Product Code	CSB-YP717293RA
Abbreviation	Ptgis
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q62969
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Source	Yeast
Target Names	Ptgis
Protein Names	Recommended name: Prostacyclin synthase EC= 5.3.99.4 Alternative name(s): Prostaglandin I2 synthase
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal His-tagged Tag-Free
Protein Length	partial
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.