



Recombinant Bacillus cereus Formamidase (amiF)

Product Code	CSB-EP717399BAAD-B
Abbreviation	amiF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q635Y7
Product Type	Recombinant Protein
Immunogen Species	Bacillus cereus (strain ZK / E33L)
Purity	≥85% (SDS-PAGE)
Sequence	MGSSGSMVKP ISGFLTALIQ YPVPVESRA DIDKQIKQII KTIHSTKAGY PGLELIVFPE YSTQQLNTHK WTTEEFLECTV PGPETDLFAE ACKESEVYGV FSIMERNPDG GEPYNTAIII DPQGEMILKY RKLNPWVPE PWKAGDLGLP VCDGPGGSKL AVCICHDGMF PEVAREAAYK GANVLIRISG YSTQVSEQWM LTNRSNAWQN LMYTLNVNLA GYDGVFYYFG EGQVCNFDGT TLVQGHRNPW EIVTAEVYPE LADQARLGWG LENNIYNLGS RGYVATPGGV KENPYTFVKD LAEGKYKVPW EDEIKVKDGT IYGYPVKKTI HS
Source	E.coli
Target Names	amiF
Protein Names	Recommended name: Formamidase EC= 3.5.1.49 Alternative name(s): Formamide amidohydrolase
Expression Region	1-332
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.