



Recombinant Bacillus cereus Methionyl-tRNA formyltransferase (fmt)

Product Code	CSB-MP717407BAAD
Abbreviation	fmt
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q636G0
Product Type	Recombinant Protein
Immunogen Species	Bacillus cereus (strain ZK / E33L)
Purity	≥85% (SDS-PAGE)
Sequence	MIKVVFMGTP DFSVPVLRRL IEDGYDVIGV VTQPDRPVGR KKVLTPTPVK VEAEKHGIPV LQPLRIREKD EYEKVLALAP DLIVTAAFGQ IVPNEILEAP KYGCINVHAS LLPRLRGGAP IHYAIMEGKE KTGITIMYMV EKLDAGDILT QVEVEIEERE TTGSLFDKLS EAGAHLLSKT VPLLIQGKLE PIKQNEEEVT FAYNIKREQE KIDWTKTGEE VYNHIRGLNP WPVAYTTLAG QVVKVWWGEK VPVTKSAEAG TIVAIEEDGF VVATGNETGV KITELQPSGK KRMSCSQFLR GTKPEIGTKL GENA
Source	Mammalian cell
Target Names	fmt
Protein Names	Recommended name: Methionyl-tRNA formyltransferase EC= 2.1.2.9
Expression Region	1-314
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.