



Recombinant Bacillus cereus Ribonuclease HII (rnhB)

Product Code	CSB-EP717408BAAD
Abbreviation	rnhB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q636I9
Product Type	Recombinant Protein
Immunogen Species	Bacillus cereus (strain ZK / E33L)
Purity	≥85% (SDS-PAGE)
Sequence	<pre>MQKVTIQEAE HLLQEIMSEE DDRFQILIKD ERKGVQKLIS KWKYKQKELAQ KEKEKFLEMS KYENALREKG LTYIAGIDEV GRGPLAGPVV TAAVILPEDF YIPGLNDSKK LSEAKRERFY DEIKAQAIAI GVGIVSPQVI DEINIYQATK QAMLDVANL SCTPEYLLID AMKLPTPIQ TSIKGDACS ISISAASIIA KVTRDRMMKE LGEKYPAYGF EQHMGYGTKQ HLEAIEVHGV LEEHRKSFAP IKDMIQK</pre>
Source	E.coli
Target Names	rnhB
Protein Names	Recommended name: Ribonuclease HII Short name= RNase HII EC= 3.1.26.4
Expression Region	1-257
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.