



Recombinant Bacillus cereus Deoxyribose-phosphate aldolase (deoC)

Product Code	CSB-EP717425BAAD
Abbreviation	deoC
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q63CR9
Product Type	Recombinant Protein
Immunogen Species	Bacillus cereus (strain ZK / E33L)
Purity	≥85% (SDS-PAGE)
Sequence	MNIAKLIDHT ILKANTTKED VMKVIEEAKE YKFASVCINP TWVKLAAEEL AGHDVDVCTV IGFPLGASTT ETKAFETKDA IAKGATEVDM VINVGALKDG DNELVEKDIY EVVQAAKGKA LVKVIETCL LTDEEKVRAC ELSVKAGADF VKTSTGFSTG GATAEDIALM RKTVGPNVGV KASGGVRTRE DAEKMVAAGA SRVGASASVA IVLNDKAGAT DNY
Source	E.coli
Target Names	deoC
Protein Names	Recommended name: Deoxyribose-phosphate aldolase Short name= DERA EC= 4.1.2.4 Alternative name(s): 2-deoxy-D-ribose 5-phosphate aldolase Phosphodeoxyriboaldolase Short name= Deoxyriboaldolase
Expression Region	1-223
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.