



Recombinant Bacillus cereus Isopentenyl-diphosphate delta-isomerase (fni)

Product Code	CSB-EP717429BAAD-B
Abbreviation	fni
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q63DN3
Product Type	Recombinant Protein
Immunogen Species	Bacillus cereus (strain ZK / E33L)
Purity	≥85% (SDS-PAGE)
Sequence	MVRAKRKLDH IEYALSTGQS RTHGFHDIDF VHQSLPNSSY ETITCETKIG ELSLSSPIFI NAMTGGGGEK TLHINEQLAY VAKHHNLAMA VGSQMAALKD ESEAASYKII RKNVPNGIFF ANLGSEATIE QAERAVDMIE ANALQIHLNV IQELTMPEGD RDFTGVLQRI EKIVLNSKVP VIVKEVGFGM SKETMQQLAS VGVTADIGG QGGTNFAAVE NERRQRMLSY FNNWGIQTAT SIEATSTNN NLSFIASGGI QTALDVAKAI ALGANTTAFAGYFLRILMED GIEKLVDIED LLHTDLKFIM TALGAKTIEE LQSVPLVIKG ETYHWLTQRG IDTTHYSRR
Source	E.coli
Target Names	fni
Protein Names	Recommended name: Isopentenyl-diphosphate delta-isomerase Short name= IPP isomerase EC= 5.3.3.2 Alternative name(s): Isopentenyl pyrophosphate isomerase
Expression Region	1-349
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.