



Recombinant *Xenopus laevis* Secretagogin (scgn)

Product Code	CSB-MP717513XBE
Abbreviation	scgn
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q63ZJ3
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	≥85% (SDS-PAGE)
Sequence	MDSAFLKLDA AEFLEIWQRF DADDNGYIEG KELDEFFCHL LKKLGITNAIT ADKVQGVKDR FMSAYDITAD GRLQIQELAN MILPEDENFL LLFRRETPLD NSVEFMRIWR KFEDEDSSGFI SAVELRNFLQ DLFLQHKKHI TGDKLDEYTD TMMKLFNRNK DGRDLNDLA KILALQENFL LQFKMDASSQ EERKSDFETI FAHYDVSKTG ALEGPEVDGF VKDMMELVKP SISGVDLDFK RQILLNHCDV NKDGKIQKSE LALCLGLKAN P
Source	Mammalian cell
Target Names	scgn
Protein Names	Recommended name: Secretagogin
Expression Region	1-271
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.