



Recombinant *Xenopus laevis* Putative deoxyribonuclease TATDN1 (tatdn1)

Product Code	CSB-BP717523XBE
Abbreviation	tatdn1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q640V9
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	≥85% (SDS-PAGE)
Sequence	MSLYRFIDIG INLTDPMFRG LYRGTRKHQD DFADIIERAV RTGVQKFMIT GGNLHESKEA IQLAQSNDRF YSTVGCHPTR CGEFEQGD PD QYLAELQNLL EDNKGKVVAV GECGLDFDRL EFC SKETQLK YFEKQFDLAE RSRLPMFLHC RNAHKEFLEI MQRNRDRCVG GVVHSFDG TK EDAEAIALD LYIGINGCSL KTESNLDVLK SIPSERLMIE TDAPWCGVKN THAGSKLVKT TFPTKKKWES GHCLKDRNEP CHIIQVLEIM ASAREEEPLE LSKTLYNNTL KLFFPAA
Source	Baculovirus
Target Names	tatdn1
Protein Names	Recommended name: Putative deoxyribonuclease TATDN1 EC= 3.1.21.-
Expression Region	1-297
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.