



Recombinant Mouse Transcription factor Dp-2 (Tfdp2)

Product Code	CSB-BP717526MO
Abbreviation	Tfdp2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q64163
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	TAKNVGLPS TNAKLRGFID QNFSPSKGNI SLVAFPVSST NSPTKILPKT LGPINVNVGP QMIISTPQRI ANSGSVLIGN PYTPAPAMVT QTHIAEAAGW VPSDRKRARE FIDSDFSESK RSKKGDKNKG GLRHFSMKVC EKVQRKGTTS YNEVADELVS EFTNSNNHLA ADSAYDQENI RRRVYDALNV LMAMNISSL PTGKKRNQVD CNSAQECQNL EIEKQRRIER IKQKRAQLQE LLLQQIAFKN LVQRNRQNEQ QNQGPPAVNS TIQLPFIIN TSRKTVIDCS ISSDKFEYLF NFDNTFEIHD DIEVLKRMGM SFGLESGKCS LEDLKIARSL VPKALEGYIT DISTGPSWLN QGLLLNSTQS VSNLDPPTGA TVPQSSVNQG LCLDAEVALA TGQLPASNSH QSSSAASHFS ESRGETPCSF NDEDEEDEEE DPSSPE
Source	Baculovirus
Target Names	Tfdp2
Protein Names	Recommended name: Transcription factor Dp-2 Alternative name(s): Dp-3 E2F dimerization partner 2
Expression Region	2-446
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.