



# Recombinant Mouse Transcription factor Dp-2 (Tfdp2)

<b>Product Code</b>	CSB-EP717526MO-B
<b>Abbreviation</b>	Tfdp2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q64163
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	TAKNVGLPS TNAKLRGFID QNFSPSKGNI SLVAFVVSST NSPTKILPKT LGPINVNVGP QMIISTPQRI ANSGSVLIGN PYTPAPAMVT QTHIAEAAGW VPSDRKRARE FIDSDFSESK RSKKGDKNKG GLRHFSMKVC EKVQRKGTTS YNEVADELVS EFTNSNNHLA ADSAYDQENI RRRVYDALNV LMAMNISSL PTGKKRNQVD CNSAQECQNL EIEKQRRIER IKQKRAQLQE LLLQQIAFKN LVQRNRQNEQ QNQGPPAVNS TIQLPFIIN TSRKTVIDCS ISSDKFEYLF NFDNTFEIHD DIEVLKRMGM SFGLESGKCS LEDLKIARSL VPKALEGYIT DISTGPSWLN QGLLLNSTQS VSNLDPPTGA TVPQSSVNQG LCLDAEVALA TGQLPASNSH QSSSAASHFS ESRGETPCSF NDEDEEDEEE DPSSPE
<b>Source</b>	E.coli
<b>Target Names</b>	Tfdp2
<b>Protein Names</b>	Recommended name: Transcription factor Dp-2 Alternative name(s): Dp-3 E2F dimerization partner 2
<b>Expression Region</b>	2-446
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.