



# Recombinant Mouse Complement component receptor 1-like protein (Cr1l), partial

<b>Product Code</b>	CSB-BP717575MO1
<b>Abbreviation</b>	Crry
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q64735
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	ELRGGLGKHG HTVHREPAVN RLCADSKRWS GLPVSAQRPF PMGHCPAPSQ LPSAKPINLT DESMFIGTY LLYECLPGYI KRQFSITCKQ DSTWTS AEDK CIRKQCKTPS DPENGLVHVH TGIQFGSRIN YTCNQGYRLI GSSSAVCVIT DQSVDWDTEA PICEWIPCEI PPGIPNGDFF SSTREDFHYG MVVTYRCNTD ARGKALFNLV GEPSLYCTSN DGEIGVWSGP PPQCIELNKC TPPPYVENAV MLSENRS LFS LRDIVEFRCH PGFIMKGASS VHCQSLNKWE PELPSCFKGV ICRLPQEMSG FQKGLGMKKE YYYGENVTLE CEDGYTLEGS SQSQCQSDGS WNPLLAKCVS RSISG
<b>Source</b>	Baculovirus
<b>Target Names</b>	Cr1l
<b>Protein Names</b>	Recommended name: Complement regulatory protein Crry Alternative name(s): Protein p65
<b>Expression Region</b>	41-405
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Partial
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.