



Recombinant Mouse Complement component receptor 1-like protein (Cr1l), partial

Product Code	CSB-MP717575MO1
Abbreviation	Crry
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q64735
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	ELRGGLGKHG HTVHREPAVN RLCADSKRWS GLPVSAQRPF PMGHCPAPSQ LPSAKPINLT DESMFIGTY LLYECLPGYI KRQFSITCKQ DSTWTS AEDK CIRKQCKTPS DPENGLVHVH TGIQFGSRIN YTCNQGYRLI GSSSAVCVIT DQSVDWDTEA PICEWIPCEI PPGIPNGDFF SSTREDFHYG MVVTYRCNTD ARGKALFNLV GEP SLYCTSN DGEIGVWSGP PPQCIELNKC TPPPYVENAV MLSEN RSLFS LR DIVEFRCH PGFIMKGASS VHCQSLNKWE PELPSCFKGV ICRLPQEMSG FQKGLGMKKE YYYGENVTLE CEDGYTLEGS SQSQCQSDGS WNPLLAKCVS RSISG
Source	Mammalian cell
Target Names	Cr1l
Protein Names	Recommended name: Complement regulatory protein Crry Alternative name(s): Protein p65
Expression Region	41-405
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.