



Recombinant *Drosophila erecta* Kinetochore protein Spc25 (Spc25)

Product Code	CSB-EP717582DLK
Abbreviation	Spc25
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q64EW1
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila erecta</i> (Fruit fly)
Purity	≥85% (SDS-PAGE)
Sequence	MAIIMAESSY ERRVKALY EK QIRMETLEGK FIKKVFKEFNS ILLDVKEAAA RHQRKVEKMQ KVVIERREEL EKRVSFMGQL AQEVEATKLR NLAMKDRIKQ QKMLARERNN EIMERIHTLS KTTGTYNQEQ ALPARVKGVT VLRGDKRDQL IPFDLNATDA EGLNSLCQHL ESLNVDVSQW QQLVSLVMDV TMEARAPTP PKEAANCKSI IEIDLTSPTS
Source	E.coli
Target Names	Spc25
Protein Names	Recommended name: Kinetochore protein Spc25
Expression Region	1-220
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.