



Recombinant Balaenoptera physalus Somatotropin (GH1)

Product Code	CSB-EP717660BED
Abbreviation	GH1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q659Q8
Product Type	Recombinant Protein
Immunogen Species	Balaenoptera physalus (Fin whale) (Balaena physalus)
Purity	>85% (SDS-PAGE)
Sequence	FPAM PLSSLFANAV LRAQHLHQLA ADTYKEFERA YIPEGQRYSI QNAQA AFCFS ETIPAPTGKD EAQQRSDVEL LRFSLLLIQS WLGPVQFLSR VFTNSLVFGT SDRVYEK LKD LEEGIQALMR ELEDGSPRAG QILKQTYDKF DTNMRSDDAL LKNYGLLS CF KKDLHKAETY LRVMKCRRFV ESSCAF
Source	E.coli
Target Names	GH1
Protein Names	Recommended name: Somatotropin Alternative name(s): Growth hormone
Expression Region	27-216
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.