



Recombinant *Bacillus licheniformis* Pantothenate synthetase (panC)

Product Code	CSB-EP717688BQU-B
Abbreviation	panC
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q65I57
Product Type	Recombinant Protein
Immunogen Species	<i>Bacillus licheniformis</i> (strain ATCC 14580 / DSM 13 / JCM 2505 / NBRC 12200 / NCIMB 9375 / NRRL NRS-1264 / Gibson 46)
Purity	≥85% (SDS-PAGE)
Sequence	MIQIYTAKDI QNLTKNYRKE GKTIGFVPTM GFLHEGHMSL VEKARKDNDI VVMSIFVNPM QFGPGEDYEA YPRDIERDRQ LAAGSGVDVL FTPEPEEMYG QEPTVTASVK KRTDVLGGRS REGHFDGVAT VLTKLFNLTS PTRVYFGMKD AQQVAVVDGL INDDFMDIEL IPVETKREED GLAKSSRVN LRAEERQEAT ALYRALQRGA ELIRNGERDP EAVKREIRSI LEKTSVIDY ADIYSYPDL IRDPLTGKVI IAVAVQFSKA RLIDNIIVDI PADQKEDL
Source	<i>E. coli</i>
Target Names	panC
Protein Names	Recommended name: Pantothenate synthetase Short name= PS EC= 6.3.2.1 Alternative name(s): Pantoate--beta-alanine ligase Pantoate-activating enzyme
Expression Region	1-288
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.