



# Recombinant *Oryza sativa* subsp. japonica

## Probable protein phosphatase 2C 51

### (Os05g0572700, LOC\_Os05g49730)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP717788OFG   |
| <b>Abbreviation</b>      | Os05g0572700, LOC_Os05g49730  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q65XK7  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | <i>Oryza sativa</i> subsp. japonica (Rice)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MRETGATDEG HACEVVVAGG DGKAAAARRR RRLELRRRLGL AAEDDAAAKR<br>IRSVKDGSSS DDSSTEVVPR SWPACVSHGS VSVIGRRREM EDAVAIERTF<br>MASTGDGAGA IRGGGEGEED FFAVYDGHGG SRVAEACRKR<br>MHVVLAAEVS LRRLRGQSAS GGDVRWKEAM LASFARM DGE<br>VVGSVAAAAP RVDGTEPSGF RTVGSTAVVA VVGRRRIVVA NCGDSRAVLS<br>RGGVALPLST DHKPDRPDEL ERVEAAGGRV INWNGYRVLG VLATRSRIGD<br>YYLKPFSVAE PEVRVVERTD KDEFLILASD GLWDVVSNEV ACKIARNCLN<br>GRAASMPES VSGSSAADAA ALLAELAVSR GSRDNISVVV VELRRLKSRA A |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | Os05g0572700  |
| <b>Protein Names</b>     | Recommended name: Probable protein phosphatase 2C 51 Short name=<br>OsPP2C51 EC= 3.1.3.16   |
| <b>Expression Region</b> | 1-381   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life   |



of lyophilized form is 12 months at -20°C/-80°C.