



Recombinant Carnation ringspot virus Movement protein (ORF3)

Product Code	CSB-EP717794CGAB-B
Abbreviation	ORF3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q66095
Product Type	Recombinant Protein
Immunogen Species	Carnation ringspot virus (isolate Lommel) (CRSV)
Purity	≥85% (SDS-PAGE)
Sequence	MAVHVDNVKD LAKTNEGVS ILNRYTDWKC KSGVTEAPLL PSNLGTRITD YAKTTLKGDS VALNFSDLVL SLGATVSVGT PGSVLVELIN PNMDGPFQLV QGQSLSWSPG SGRPCLMIFS IHHQLTADAE PFRIRISNNG IPTKKTFCARC HAYWGFDLSP RMRYYKNEPA KRIDLVDVGFY KTHLSNMKQV RDYVQYTFDN SRMDGNPQLV AKSTMNVVPR IADVPKYVGI APPSRSGNHQ EATPDDWLKQ YVDKDSETNK LSDVESSSDS SLLSMRARS RRYTKNYKIP IKRHPQATGE VGTT
Source	E.coli
Target Names	ORF3
Protein Names	Recommended name: Movement protein Short name= MP
Expression Region	1-304
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.