



Recombinant Protein FdhD homolog (fdhD)

Product Code	CSB-BP717824YAH
Abbreviation	fdhD
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q663U8
Product Type	Recombinant Protein
Immunogen Species	Yersinia pseudotuberculosis serotype I (strain IP32953)
Purity	≥85% (SDS-PAGE)
Sequence	MSQIKPSRLS SSAEIRGARQ LDVLQRHKLA EPQQDWLAEE VPVALVYNGI SHVMMATPK DLAAFALGFS LSEGISSPQ DIYAIEMTPG CNGIEVNIEL SSRRFAGLKE RRRAMAGRTG CGVCGIEQLD DIFRPITPLP FTQAFNLEHL DTALAQLKQV QPVGQLTGCT HAAAWINPEG ELLGGCEDVG RHVALDKLLG IRAKQPWQQG AVLVSSRASY EMVQKTAMCG AEILFAVSAA TTLAVEVAER CNLTLVGFSK PGRATVYTHP QRIK
Source	Baculovirus
Target Names	fdhD
Protein Names	Recommended name: Protein FdhD homolog
Expression Region	1-274
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>