



# Recombinant DNA-directed RNA polymerase subunit alpha (rpoA)

<b>Product Code</b>	CSB-EP717834YAH-B
<b>Abbreviation</b>	rpoA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q664U6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Yersinia pseudotuberculosis serotype I (strain IP32953)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MQGSVTEFLK PRLVDIEQVS STHAKVTLEP LERGFHTLG NALRRILLSS MPGCAVTEVE IDGVLHEYST KEGVQEDILE ILLNLKGLAV RVQ GKDEVIL TLNKSGIGPV TAADITHDGD VEIVKPQHVI CHLTDENASI NMRIKVQRGR GYVPASARIH SEEDERPIGR LLVDACYSPV ERIAYNVEAA RVEQRTDLDK LVIEMETNGT IDPEEAIIRA ATILAEQLEA FVDLRDVRQP EVKEEKPEFD PILLRPVDDL ELTVRSANCL KAEAIHYIGD LVQRTEVELL KTPNLGKKS L TEIKDVLASR GLSLGMRLLEN WPPASIADE
<b>Source</b>	E.coli
<b>Target Names</b>	rpoA
<b>Protein Names</b>	Recommended name: DNA-directed RNA polymerase subunit alpha Short name= RNAP subunit alpha EC= 2.7.7.6 Alternative name(s): RNA polymerase subunit alpha Transcriptase subunit alpha
<b>Expression Region</b>	1-329
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.