



Recombinant Equine herpesvirus 2 Gene 55 protein (55)

Product Code	CSB-BP717850EFP
Abbreviation	55
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q66657
Product Type	Recombinant Protein
Immunogen Species	Equine herpesvirus 2 (strain 86/87) (EHV-2)
Purity	≥85% (SDS-PAGE)
Sequence	MSSRWIQWTC CGLWPFGKSS RQYRQIINVY EDDSERVRAE IDMG LPPGVS VGDLLQNRQG EEAIRQAHML ALQSNAITDY LQRFDSVRVP ESCQSIVSAQ IAKLRSVQHI IWNTMISMAV GNV TIDDSTL KVLLNKQAGE SMALMEMEKL AAALKMDGST DWAQNIASVV TAPPRAAREQ MAFREPIYDD PEGGMEMIDL SSTPGPTPAP TPAPSMERVA VPAS
Source	Baculovirus
Target Names	55
Protein Names	Recommended name: Gene 55 protein
Expression Region	1-224
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.