



Recombinant Glucokinase (glk)

Product Code	CSB-EP717867YAH-B
Abbreviation	glk
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q668N7
Product Type	Recombinant Protein
Immunogen Species	Yersinia pseudotuberculosis serotype I (strain IP32953)
Purity	≥85% (SDS-PAGE)
Sequence	MTTYALVGDV GGTNARLALC AVATGEILQA KTYSGLEYES LEDVIKQYLS EHQAKVTDAC IAIACPITGD WVAMTNHTWA FSIAAMQQNL GLDHLEVIND FTAVSMAIPV LPAQDVLQFG GTQPQPGKPV AVYGAGTGLG VAHLVNVDRR WISLAGEGGH VDFAPNSEEE DQILAVLRQE LGHVSAERVL SGPGLVNLVYR AIVISDARLP EKLAPKDITA RALADSCTDC RRALSLFCVI MGRFGGNLAL NLSTFGGVYI AGGIVPRFME FFKASGFRAA FEDKGRFKDF LQDIPVYMIT HPQPGLLGAG AYLRQKLGYE LSS
Source	E.coli
Target Names	glk
Protein Names	Recommended name: Glucokinase EC= 2.7.1.2 Alternative name(s): Glucose kinase
Expression Region	1-323
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>