



Recombinant Acetate kinase (ackA)

Product Code	CSB-EP717869YAH
Abbreviation	ackA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q668Z0
Product Type	Recombinant Protein
Immunogen Species	Yersinia pseudotuberculosis serotype I (strain IP32953)
Purity	≥85% (SDS-PAGE)
Sequence	MSSKLVLVLN CGSSSLKFAI IDATNGEEHI SGLAECFHLP EARIKWKVGD GKQEAALGAG AAHSEALNFI VNTILAQKPA LSAQLTAIGH RIVHGGEKFT SSVIVTDDVI QGIKDSIPFA PLHNPAHLIG IAEALKSFPN LADKNVAVFD TAFHQTMPEE SYLYALPYSL YKDHGIRRYG AHGTSHFYVS QEAAILNKP LEELNVITCH LGNGGSVTAV RNGKCVDTSM GLTPLEGLVM GTRSGDLDP A IIFHLHDAMG MSVDQINTLL TKESGLLGLT EVTSDCRYVE DNYATKADAK RAMDVFCHRL AKYIGSYTAL MDGRLDAVVF TGGIGENAAM VRELTLDKLG LLGFEIDHER NMAARFGKSG TITKDSSRLA LVIPTNEELV IAQDAARLTA
Source	E.coli
Target Names	ackA
Protein Names	Recommended name: Acetate kinase EC= 2.7.2.1 Alternative name(s): Acetokinase
Expression Region	1-400
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.