



# Recombinant Ribosomal RNA large subunit methyltransferase F (rlmF)

<b>Product Code</b>	CSB-EP717872YAH-B
<b>Abbreviation</b>	rlmF
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q669D4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Yersinia pseudotuberculosis serotype I (strain IP32953)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	<p>           MLSYAPENAY QRASTMENKK VFPKEKSLGH PRNRHRSRYD FDALSVSCPE            LIPFLAPTAY GDISVDFADP LAVKMLNKAL LKHFYGIIEW DIPADSLCPP            IPGRADYVHH LADLLASCNG EVIPKGKNIA LLDIGVGANC IYPIIGQREY            GWRFTGTDID SHALSAKMV VSMNPTLKNT LRLKQQKDPH AIFEGVWAVN            ERYDATLCNP PFHGSAAEAA ATTRRKLHKL GKNEVAAPV QNFGGKNSL            WCEGGEFV SRMVAESVAK AQNCFWFTSL ISKKTLPAL YHALRYVKAV            EVRTIEMAQG QKVSFRVAWT FLTPEQQAAW VAERWA         </p>
<b>Source</b>	E.coli
<b>Target Names</b>	rlmF
<b>Protein Names</b>	Recommended name: Ribosomal RNA large subunit methyltransferase F EC=2.1.1.181 Alternative name(s): 23S rRNA mA1618 methyltransferase rRNA adenine N-6-methyltransferase
<b>Expression Region</b>	1-336
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.