



Recombinant Imidazolonepropionase (hutI)

Product Code	CSB-EP717883YAH
Abbreviation	hutI
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q66B14
Product Type	Recombinant Protein
Immunogen Species	Yersinia pseudotuberculosis serotype I (strain IP32953)
Purity	≥85% (SDS-PAGE)
Sequence	MVSVTHCDSL WFGADIITMR GGN YQLIPQG AIAVTGDKIV WIGPHAE LPP IHAARQVVYE GGLITPGLID CHTHLVFGGD RSNEFEQRLN GVS YAEIAAN GGGIISTVRA TRQASEQQLL EQALFRLKPL LAEGVT TIEI KSGYGLNLES EIKMLRVARR LGELLPIDVK TTCLAAHALP PEFIGQPDDY IDVVCNSIIP QVAVENLADA VDAFCEHLAF SPAQVERVFL AAQKAGLPVK LHAEQLSALR GATLAAKFHA ISADHLEYAT ESDVQAMAKA GTVAVLLPGA YYLLRETQCP PIDLFRQYKV PMALASDANP GTSPVLSLRL MLNMACTLFR MTP EEALAGV TCHAAQALGV QQTQGTLETG KLANWVHWPL SHPAELAYWL GGQLPATVVF RGEVRP
Source	E.coli
Target Names	hutI
Protein Names	Recommended name: Imidazolonepropionase EC= 3.5.2.7 Alternative name(s): Imidazolone-5-propionate hydrolase
Expression Region	1-406
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.