



# Recombinant Imidazolonepropionase (hutI)

<b>Product Code</b>	CSB-YP717883YAH
<b>Abbreviation</b>	hutI
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q66B14
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Yersinia pseudotuberculosis serotype I (strain IP32953)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MVSVTHCDL WFGADIITMR GGNYQLIPQG AIAVTGDKIV WIGPHAE LPP IHAARQVVYE GGLITPGLID CHTHLVFGGD RSNEFEQRLN GVSYAEIAAN GGGIISTVRA TRQASEQQLL EQALFRLKPL LAEGVTTIEI KSGYGLNLES EIKMLRVARR LGELLPIDVK TTCLAAHALP PEFIGQPDDY IDVVCNSIIP QVAVENLADA VDAFCEHLAF SPAQVERVFL AAQKAGLPVK LHAEQLSALR GATLAAKFHA ISADHLEYAT ESDVQAMA KA GTVAVLLPGA YLLRETQCP PIDLFRQYKV PMALASDANP GTSPVLSRL MLNMACTLFR MTPPEALAGV TCHAAQALGV QQTQGTLETG KLANWVHWPL SHPAELAYWL GGQLPATVVF RGEVRP
<b>Source</b>	Yeast
<b>Target Names</b>	hutI
<b>Protein Names</b>	Recommended name: Imidazolonepropionase EC= 3.5.2.7 Alternative name(s): Imidazolone-5-propionate hydrolase
<b>Expression Region</b>	1-406
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.