



Recombinant Deferrochelataase/peroxidase EfeB (efeB)

| | |
|--------------------------|--|
| Product Code | CSB-MP717884YAH |
| Abbreviation | efeB |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q66BP4 |
| Product Type | Recombinant Protein |
| Immunogen Species | Yersinia pseudotuberculosis serotype I (strain IP32953) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | QAADC RSPDVAGTQD ERWQKQPFYG QHQAGVLTPQ QAAMMLVAFD VLATDKTSLI RLFKLLTERL AFLTTGGRAP SVNAKLPPLD SGIMGPEIYP DNLTVTVSVG NALFDERFGL QGQKPLRLQR MTRFPNDSL D AGLCHGDVML QICANTNETV IHALRDIKH TPDLLSVRWK REGFISAHAA RSKGQETPIN LLGFKDGTAN PKISNKPLIN NVVWVSNNAG EPAWAVGGSY QVVRIIRFKV EFWDRTPHQE QQTIFGRDKN SGAPLGMQHE HDEPNYAKDP EGKVIPMDAH IRLANPRTIE TQRNMLRRG YSYSLGVSNS GQLDMGLLFV CYQSDLAQAF LTVQERLNGE ALEEYVKPIG GGYFFTLPGV ADANHILAQS LLEA |
| Source | Mammalian cell |
| Target Names | efeB |
| Protein Names | Recommended name: Deferrochelataase/peroxidase EfeB EC= 1.11.1.- |
| Expression Region | 46-434 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |