



Recombinant Thiamine-phosphate synthase (thiE)

Product Code	CSB-YP717926YAH
Abbreviation	thiE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q66FP6
Product Type	Recombinant Protein
Immunogen Species	Yersinia pseudotuberculosis serotype I (strain IP32953)
Purity	≥85% (SDS-PAGE)
Sequence	MATPGFPSTE QRLGLYPVVD SLLWIERLLA AGVTTLQLRI KNADDAQVEQ DIVAAIELGK RYQARLFIND YWQLAVKHGA YGVHLGQEDL ETADLAAIQQ AGLRLGISTH DEHELAVAKT LRPSYIALGH IFPTQTKQMP SSPQGLASLS RQVKNTPDYP TVAIGGISIE RVPHVLATGV GSVAVVSAIT LASDWQRATA QLLHLIEGKE LADEK
Source	Yeast
Target Names	thiE
Protein Names	Recommended name: Thiamine-phosphate synthase Short name= TP synthase Short name= TPS EC= 2.5.1.3 Alternative name(s): Thiamine-phosphate pyrophosphorylase Short name= TMP pyrophosphorylase Short name= TMP- PPase
Expression Region	1-215
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.