



Recombinant ATP-dependent protease subunit HslV (hslV)

Product Code	CSB-BP717933YAH
Abbreviation	hslV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q66G84
Product Type	Recombinant Protein
Immunogen Species	Yersinia pseudotuberculosis serotype I (strain IP32953)
Purity	≥85% (SDS-PAGE)
Sequence	TTIVSVRRD GHVVIGGDGQ VTLGNTVMKG NAKKVRRLYN NKVIAGFAGG TADAFTLFEL FERKLEMHQG HLTKAAVELA KDWRTDRMLR KLEALLAVAD ETASLIITGN GDVVQPEDDL IAIGSGGPYA QSAARALLEN TELGARDIVE KLSIAGDIC IYTNRFQTIE ELTY
Source	Baculovirus
Target Names	hslV
Protein Names	Recommended name: ATP-dependent protease subunit HslV EC= 3.4.25.2
Expression Region	2-174
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.