



# Recombinant Rat Protein FAM71B (Fam71b)

<b>Product Code</b>	CSB-BP717939RA
<b>Abbreviation</b>	Fam71b
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q66H38
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSSECLLPYY TAHSYRSMGV FNTSMGNLQR QLYKGGEYDI FKYAPMFESD FIQISKKGEV IDVHNRVRMV TVCIASSTSPV LPLPDVMLLA RPAKVCEEHA RRARFIKGRG YKPSKTLELT RLLPLKFVKI SIHDHEKQQL RLKLATGRTF YLQLCPSSDA REDLFCYWEEK LVYLLRPPVN SCISNPSIPT ADTSTETKST LVSEIHGEGD RDSKFQTSQD VSEATSAAFA GGERTQPANIA AILTTGLAKA RAAGAAAGTG AAGTAGTAGA AGATGAAGAT GAAGSARAAG GAGSARAAGA TGSARAAGAT GSARAAGATG SARAAGGNTA AAVMTPLAGL ARAGTPTSPV SGVISIAATT SKSPGSGQVA TGLTGTASKD QERSESSKAM AVVANITTES VDVVLAGAAS FTSESPSAEG DAYGSPDTGL NVAFAGSITT KGPAEDKPEA PLVSTLQSEG YMCERDGSQK VSHTSSETQK EKRRERRESDR KGSRKSSHHQ RTGASRHSSS KDKGRKTSSY RSVSGHGKTR EDKKEKGRGS LRDQRHSSSY RSESRTGHKS RKNRPAASGG FVSKRATKIR SFFRAFLVRP TLKTENTSGE RGGVDIVTKL VEKQDIETTV EKSKDLEFSD TMASESMEKI ILETKPI
<b>Source</b>	Baculovirus
<b>Target Names</b>	Fam71b
<b>Protein Names</b>	Recommended name: Protein FAM71B
<b>Expression Region</b>	1-647
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.