



Recombinant Hepatitis B virus genotype C subtype adr Putative X-Core fused protein

Product Code	CSB-MP717989HED
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q67863
Product Type	Recombinant Protein
Immunogen Species	Hepatitis B virus genotype C subtype adr (isolate Korea/Kim/1989) (HBV-C)
Purity	≥85% (SDS-PAGE)
Sequence	MAARVCCQLD PARDVLCLRP VGAESRGRPV SRPFGPHPS P SSSAVPADHG AHLRLGLPV CAFSSAGPCA LRFTSARRME TTVNADQVLP KVLHKRTLGL SAMSTSDLEA YFKDCLFKDW EDLGEEIRLM IFVLGGCRHK LVCSPAPCNF FHLCLIISCS CPTVHASKLC LGWLWGMDID PYKEFGASVE LLSFLPSDF PSIRDLLDTA SALYREALS PEHCSPHHTA LRQAILCWGE LMNLATWVGS NLEDPASREL VVSYYNVNMG LKIRQLLWFH ISCLTFGRET VLEYLVSGFV WIRTPPAYRP PNAPILSTLP ETTVVRRRGR SPRRRTPSPR RRRSESPRRR RSQSRESQC
Source	Mammalian cell
Protein Names	Recommended name: Putative X-Core fused protein
Expression Region	1-359
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.