



# Recombinant Acetate kinase (ackA)

<b>Product Code</b>	CSB-BP718024SRW
<b>Abbreviation</b>	ackA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q67L76
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Symbiobacterium thermophilum (strain T / IAM 14863)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MKILVLNCGS SSIKYQLLDM ETESVLAIGK VERVGMEDAI FEYKGAGEKI KEIGTILDHT AGVRKVLKAL TDPVHGVIIQA PEEIAAVGHR VVHGGEAFVK SVIITDEVKR TLRELFDLAP LHNPNALGTI LAEEAALPHA PGVAVFDTAF HATMPRKAFL YALPYAMYKG HKVRRYGFHG TSHKYVSQRA ADLLGRPLEE VNLITCHLGN GSSVTAVQGG KSVDTSMGFT PLAGLIMGTR AGDLDPGVIP YIMAREEIGI SEVNSMLNKH SGILGISGIS SDMRQIEEEM EKGDPRAVET FEMMEYRLRK YIGAYAALVG RVDGIVFTGG IGENSPLFRK AVCQNLGFLG VTFDEEANNC RGQERIISGP DSKVAVMVIP TNEELMIARE TAEIIRNR
<b>Source</b>	Baculovirus
<b>Target Names</b>	ackA
<b>Protein Names</b>	Recommended name: Acetate kinase EC= 2.7.2.1 Alternative name(s): Acetokinase
<b>Expression Region</b>	1-398
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.