



Recombinant UPF0271 protein STH2522 (STH2522)

Product Code	CSB-MP718025SRW
Abbreviation	STH2522
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q67LD9
Product Type	Recombinant Protein
Immunogen Species	Symbiobacterium thermophilum (strain T / IAM 14863)
Purity	≥85% (SDS-PAGE)
Sequence	MGAIRGPVRK IDLNCDMGES FGVYRIGADE EIMPLITSAN IACGFHGGDP QVMRRTVRLA REHGVAVGAH PGYRDLVGFWR RRPVRCSPDE VYADVLYQIG ALAAFRCRAEG VALRHVKPHG ALYNTAAADA AIAGAVARAV ADFDRSLMLY APPGSALEQA GLAAGLRVIR EGFADRGYAA DGTLLPRTHP GAVLHEPERA AAQARGMVCS GTVTADTGEE VAVPAETLCV HGDHPSVIQV LRRIRSELEA AGVSVGAPGA D
Source	Mammalian cell
Target Names	STH2522
Protein Names	Recommended name: UPF0271 protein STH2522
Expression Region	1-261
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.