



Recombinant Bifunctional protein FoD (foD)

Product Code	CSB-BP718029SRW
Abbreviation	foD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q67NB1
Product Type	Recombinant Protein
Immunogen Species	Symbiobacterium thermophilum (strain T / IAM 14863)
Purity	≥85% (SDS-PAGE)
Sequence	MTAGVIDGKQ VAAQVQEEVA REVAALKAEG LTPGLAVVLV GDDPASKVYT ANKERTAKEL GMHSVLYHLP ATTTQQELAD LVMQLNRDPA IHGILVQSPL PEGLDMDAVI RLIDPRKDVD GFHPENVGRL WIGQDGLVPC TPAGVMRLLD AYGVDPKGKH AVVVGRSNIV GKPMALLLQ RHATVTICHS RTPDLAETCR RADILVA AVG RLQMITAEHV KPGAVVIDVG INPVPGFKKR IRGDVDFDSV QEVAGLITPV PGGVGPMPTIA MLMANTVKAA RMQTGR
Source	Baculovirus
Target Names	foD
Protein Names	Recommended name: Bifunctional protein FoD Including the following 2 domains: Methylenetetrahydrofolate dehydrogenase EC= 1.5.1.5 Methenyltetrahydrofolate cyclohydrolase EC= 3.5.4.9
Expression Region	1-286
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.