



Recombinant Proline--tRNA ligase (proS)

Product Code	CSB-BP718035SRW
Abbreviation	proS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q67PA6
Product Type	Recombinant Protein
Immunogen Species	Symbiobacterium thermophilum (strain T / IAM 14863)
Purity	≥85% (SDS-PAGE)
Sequence	MKEEKSFVKE ITPQSEDFSR WYIDVIRKAD LMDYTPVRGC IVFKPDGYEL WERIQAGLDK RFKETGHRNA YFPM LIPESF FQKEKEHVEG FNP ELPWVTE AGGEKLEERL ALRPTSETII GHMYAQW IQS YRDL PVLINQ WANVFRWEKR TLPFLRTSEF LWQEGHTAHA TEEEAREETM RMLEVYRDFV ETEMAIPLYV GQKTPSEKFA GAVDTYSIEA MMKDGKALQA GTSHYLGQNF AKGFEIKFLD RDNQLKYVHT TSWGVSTR LI GALIMVHGDD RGLALPPRLA PIQVIMIPVG PPKFRDKVMA RFDPLFDALK AAGVRVKADL REETPGWKFN EWEMRGVPLR IEIGPRDVDN KQCLMVRRTD GEKIPVPLDE AVERVQALLE EIQRNMFVKA KEFRDAHSHL HINTLAELKA HIARCEETGE IAGFVLAGWC GDAACEAKVK EETKFTSRNI PFHPPARKQV CLCCGKEAQH TVWFARAY
Source	Baculovirus
Target Names	proS
Protein Names	Recommended name: Proline--tRNA ligase EC= 6.1.1.15 Alternative name(s): Prolyl-tRNA synthetase Short name= ProRS
Expression Region	1-488
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.