



Recombinant Ribonuclease HII 1 (rnhB1)

Product Code	CSB-EP718036SRW
Abbreviation	rnhB1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q67PD4
Product Type	Recombinant Protein
Immunogen Species	Symbiobacterium thermophilum (strain T / IAM 14863)
Purity	≥85% (SDS-PAGE)
Sequence	MKRLQALMAE QGPDIIPEVI AALKDDPREG AQRLVRQCQT YLKERERQRA VLQRMYNRYER QLWSMGYRHI AGIDEVGRGP LAGPVVAAAV ILPGLEVELPG IEEAKRLSDR RRMELYRRIR EVAVAVGVGL VQPEGIDEAS VVVATYKAMM KAVADLPVTP DYLLIDGVHL PNVSQPQVPV VGGETLSCSI AAAAVVAKVV RDEYMVEMDA KYPAYGFAHH KGYGTPEHRL ALERYGPCPI HRKPSDGTPA GAVLLFTEG
Source	E.coli
Target Names	rnhB1
Protein Names	Recommended name: Ribonuclease HII 1 Short name= RNase HII 1 EC= 3.1.26.4
Expression Region	1-259
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.