



# Recombinant *Oryza sativa* subsp. japonica Putative B3 domain-containing protein Os06g0632500 (Os06g0632500, LOC\_Os06g42630)

<b>Product Code</b>	CSB-BP718056OFG
<b>Abbreviation</b>	Os06g0632500, LOC_Os06g42630
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q67VL7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. japonica (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MPLSAGSLRR RPKFAAPLLS PACLHRLSVP GEFAARLDDD DAAGVGEEEEE EESGRRAAAV LVVGPLGKVV RVELRRSPAG DGEAWLGGGW SELAAAHGLG EGWGVVLRLE RRGVASLRVF DPGFCLARFC TPHAGMRTKD RPRFIKLLQQ EDLEKMKIPE KVVQQHLTET YTNNHQNAII VCPLGKFWRV ELQREQPDVL LRDGWAPFLA AHDLSEGNIL LFRYEGNMVF TVEVFLQNGC LKEYKTAALY LTDGTEGPSN APQQSAAKVG VSPVKRKRTR RIEGTCLEGP NRKSRASPIS VKVEPHKKHV SIVSQNSFTK EMTAYSIHSL LSVRGTFCSQ IGLLEACAIT LKISMKKKGS WRVAFKTANT YGYINGPGWR KFCLENELEV QRR
<b>Source</b>	Baculovirus
<b>Target Names</b>	Os06g0632500
<b>Protein Names</b>	Recommended name: Putative B3 domain-containing protein Os06g0632500
<b>Expression Region</b>	1-393
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.