



Recombinant *Oryza sativa* subsp. japonica Tubby-like F-box protein 10 (TULP10)

Product Code	CSB-YP718069OFG
Abbreviation	TULP10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q688Y7
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MAAVREPREE AAVGEGEGEE EGRWGGLLPE LVEEVVRRVE ASGGERWPAR KDLVSCACVC RRWREAAAV VRPLPESGRI TFPSSLKQPG PKDFPIQCFV KRNKKKSMFY LYLGLLNGTM DKGKFLMAAR RFRRGPHTHEY VISLDADDLS QGSNAYVGKL RSDFWGTNFK IYDNQPPYDD AKTSSTRSSQ RFGSTHFRFGS RRICPQISAG NFNVGQISYK YNLLKSRGPR RMFCTMECPS TQETWENSLK TKSLRCTGTT VLRNKAPRWH EHLQCWCLNF HGRVTVASVK NFQLVATADP SHPDSVGDEE TVILQFGKVD SNIFTMDYRQ PLSAFQAFAI CLSSFGTKLA CE
Source	Yeast
Target Names	TULP10
Protein Names	Recommended name: Tubby-like F-box protein 10 Short name= OsTLP10 Alternative name(s): Tubby-like F-box protein 11 Short name= OsTLP11
Expression Region	1-352
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.