



Recombinant Human GTPase SLIP-GC (C8orf80)

Product Code	CSB-EP718072HU-B
Abbreviation	C8orf80
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q68CJ6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAETKDVFQQ EPHPVEDDLY KERTRKRRKS DRDQRFRAFQ SMEQSALKEY EKLESRTRRV LSNTYQKLIQ SVFLDDSDIPN GVKYLINRLL ALIEKPTVDP IYIALFGSTG AGKSSLINAI IQQAMFLPVS GESICTSCIV QVSSGCCVQY EAKIHLLSDQ EWREELKNLT KLLHRTEELS REEADAWNDR EAVEEATWKL QMIYGNGAES KNYEELLRAK PKRKIPTSRV ITLKAEAEAE LSIKLDPYIR TQRRDWDGEA AEMRIWPLIK HVEVTLPKSD LIPEGVVLVD IPGTGDFNSK RDEMWWKTID KCSVIWVISD IERVSGGQAH EDLLNESIKA CQRGFCRDVA LVVTKMDKLH LPEYLRERKA GNQAIQSQRV AVLERNEMIK LQRTRILKEK LKRKLPADFK VLEASDLVYT VSAQEYWQQA LLTEEETEIP KLREYIRKSL LDKKKRTVTK YVTEAFGLLL LTDSFNSTQN LPNEHLHMSV LRRFAEEKVE LLEKAIQCF ACMEQPLQEG VRTARTSYRC ILRACLVRSK GNQGFHQTLK AVCLKNGIYA SRTLARIDLN EALTQPVYDQ IDPVFGSIFR TGKPTGSALM PHIDAFKQSL QEKMTIEGIR SGWKYDSCKK NFLIQEISAI LGGLEDHILR RKRRRIESLT ASVQSDLKLC YEEAAQITGK KACERMKDAI RRGVDRQVAE GMFERAQERM QHQFQQLKTG IVEKVKGSIT TMLALASSQG DGLYKELADV GSEYKEMEKL HRSLREVAEN ARLRKGMQEF LLRASPSKAG PPGTSL
Source	E.coli
Target Names	NUGGC
Protein Names	Recommended name: GTPase SLIP-GC EC= 3.6.1.- Alternative name(s): Speckled-like pattern in the germinal center
Expression Region	1-796
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.