



Recombinant Human GTPase SLIP-GC (C8orf80)

| | |
|--------------------------|---|
| Product Code | CSB-EP718072HU-B |
| Abbreviation | C8orf80 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q68CJ6 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MAETKDVFGQ EPHPVEDDLY KERTRKRRKS DRDQRFRAFSP SMEQSALKEY EKLESRTRRV LSNTYQKLIQ SVFLDDSSIPN GVKYLINRLL ALIEKPTVDP IYIALFGSTG AGKSSLINAI IQQAMFLPVS GESICTSCIV QVSSGCCVQY EAKIHLLSDQ EWREELKNLT KLLHRTEELS REEADAWNDR EAVEEATWKL QMIYGNGAES KNYEELLRAK PKRKIPTSRV ITLKAEAEAE LSIKLDPYIR TQRRDWDGEA AEMRIWPLIK HVEVTLPKSD LIPEGVVLVD IPGTGDFNSK RDEMWKKTID KCSVIWVISD IERVSGGQAH EDLLNESIKA CQRGFCRDVA LVVTKMDKLH LPEYLREKKA GNQAIQSQRE AVLERNEMIK LQRTRILKEK LKRKLPADFK VLEASDLVYT VSAQEYWQQA LLTEETEIP KLREYIRKSL LDKKKRTVTK YVTEAFGLLL LTDSFNSTQN LPNEHLHMSV LRRFAEEKVE LLEKAIQCF ACMEQPLQEG VRTARTSYRC ILRACLVRSK GNQGFHQTLK AVCLKNGIYA SRTLARIDLN EALTQPVYDQ IDPVFGSIFR TGKPTGSALM PHIDAFKQSL QEKMTIEGIR SGWKYDSCCK NFLIQEISAI LGGLEDHILR RKRRRIESLT ASVQSDLKLC YEEAAQITGK KACERMKDAI RRGVDRQVAE GMFERAQERM QHQFQQLKTG IVEKVKGSIT TMLALASSQG DGLYKELADV GSEYKEMEKL HRSREVAEN ARLRKGMQEF LLRASPSKAG PPGTSL |
| Source | E.coli |
| Target Names | NUGGC |
| Protein Names | Recommended name: GTPase SLIP-GC EC= 3.6.1.- Alternative name(s): Speckled-like pattern in the germinal center |
| Expression Region | 1-796 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final |



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.