



Recombinant Human U3 small nucleolar RNA-associated protein 25 homolog (UTP25)

Product Code	CSB-BP718073HU
Abbreviation	DIEXF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q68CQ4
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGKRGSRSSQS QLLNLTLLKKQ KKHLRDFGEE HPFYDRVSRK EAKPQICQLS ESSDSSDSES DSESEPQQVS GYHRLATLK NVSEEEEEDE EEEEEEDSIV DDAEMNDEDG GSDVSVEEEM AAESTESPEN VALSADPEGK EDGEPPGTS QTSPEEFTDA KHESLFSLET NFLEEEESGDN SSLKASQDPF LQHVNKELKE KAIQAVATNP KTTHELKWPI LGQLFFSSKF QKLETFKPPK DIDLKSLHLQ KPLESTWTKT NSQFLSGPQK SSSPFTPLQK ELFLIMNSYR DLFPERTAL KNGEEIRHVY CLHVINHILK ANAQVLGNNS RRRSQKFGVG DDDDFRDQGL TRPKVLIVVP FREAAALRVVQ LFISLLEGDS KKKIIVSNKK RFQGEYGSDP EERPPNLKRP EDYEAUVFVGN IDDHFRIGVA ILQRSIRLYA PFYSSDILIA SPLGLRTIIG GEGEKKRDFD FLSSIELLII DQADIYLMQN WEHVLHLMNH MNLLPLDSHG VDFSRVRMWS LNNWSKYRQ TLLFGALQDA QINSVFNKYC VNMQQQVAVR NVPMTGSISH VLVQLPHVFQ RMEAENLASV IDARFNFFVN KILPQYRDAV MSHTLIYIPS YFDFVRLRNY FKKEELNFTH ICEYTQKSGV SRARHFFLQG EKQFLLFTER FHFYKRYTIK GIRNLIFYEL PTYPHFYSEI CNMLRATNRG EEATWTCTVL YSKYDAQRLA AVVGVERAAQ MLQSNKNVHL FITGEK
Source	Baculovirus
Target Names	DIEXF
Protein Names	Recommended name: Digestive organ expansion factor homolog
Expression Region	1-756
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.