



# Recombinant *Xenopus laevis* Alpha-ketoglutarate-dependent dioxygenase FTO (fto), partial

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-BP718083XBE                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Abbreviation</b>      | fto                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                                                                                                    |
| <b>Uniprot No.</b>       | Q68F54                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Immunogen Species</b> | <i>Xenopus laevis</i> (African clawed frog)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Sequence</b>          | MKKRGVLC EE EKEVKKRKLL EQIGDGHLPY LSSKDDTFYD LWGSCYSKLT<br>LLQAKHVPVD LHHTVQNAFL SLLRNGCLFQ DLVRLKGKDI ITPVSRILIG<br>RPGYTYKYLN TRLFAVPWIN DELNTQYSTR NLLETYKA FS DLNDFLYSQT<br>VNELQKLGKI HKDNFFQHSE YKEQIKSDQR PSTSNYENVK LHSMESFNVT<br>LINYMDPQNM SFLKEEPYFG MGKMAVSWHH DENLVEQSTV AVYNYSYQES<br>SNDVNGEDED PTRWHVGLKI AWDIETPGLL MPLSSGDCYL MLDDL NKTHQ<br>HCVIAGCQPR FSSTHRVAES STGTLQYIKS QCNSALQNLH MNPDTGAADL<br>KNLEPNVLGQ TEEIHNEVEF EWLRQFWFQG KRYNKCTTFW KEAMTELEIH<br>WKQMETMTSL VLKAIENENL TVDEKCNILK NILPPLVERQ DLRHDWRERC<br>RSKLAKMLPP DQVPCFYPIW NDDNKTMLPL FDLHNIISAL QKTLEENELK L |
| <b>Source</b>            | Baculovirus                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Target Names</b>      | fto                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Protein Names</b>     | Recommended name: Alpha-ketoglutarate-dependent dioxygenase FTO EC=1.14.11.-                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Expression Region</b> | 1-501                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Protein Length</b>    | Partial                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                                                                                                                                           |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life                                                                                                                                                                                                                                                                                                                                                                     |



of lyophilized form is 12 months at -20°C/-80°C.