



Recombinant Mouse Protein SPT2 homolog (Spty2d1)

Product Code	CSB-BP718087MO
Abbreviation	Spty2d1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q68FG3
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MDFREILLIA SKGQGVNHVP KRYSLAVGPP KKDPKVKGVQ SAAVQAFLRR KEEELRQKAL EEKKRKEELV KKRIELKHDK KARAMAKRTK DNFHGYDGIP VEEKTKKKQL VESHLNQGTQ QEYDVEEEDF IDYNQAELDQ DYEEEQEPPK AESKPKAPLK SAPSPMNFTD LLRLAEKKQF EPVEIKVVKK AEDRPLTAE LREREFLEK HRKKKPEPDA KLPPPVKKKA PSHKDIMGTK PSRGAGDRQL ASKGLPFPQA EKKFRPSTAS EKQAALSSPK SLPGERTKVG SGSSTQPSLR EGHNRPVFNG AGKPRPSTCS PSVPKTPASG TQKSASEHKA KKPLPSHPSH SKPGPTVLSH NKS KSPGVRQ PGSNSGSAPG QPNPGTARPT LSSGPVPRRQ NGSSSSGPEQ SAGGIRKLAS NSHLSGRTLN GTNGPGRPAS SSSGPGRPIS GSAGSGRPVG SSGGPGQPVN NPHDLRRPMN SLSSPGRAVS GPGRSISGSI PAGRTVNSGP GRPVSSLGPG RAVSNPGLPT KPRCTVVSET ISSKNIISRS SNGQINGMKP LLSGYRSAQG PQRLPFPTGY KRPREYEEDD DDEYDSEMDD FIEDEGEPQE EISKHIREIF GYDRKKYKDE SDYALRYMES SWKEQQKEEA KSLRLGMQED LEEMRREEEE LKRRKAKKLK RH
Source	Baculovirus
Target Names	Spty2d1
Protein Names	Recommended name: Protein SPT2 homolog Alternative name(s): SPT2 domain-containing protein 1
Expression Region	1-682
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.