



Recombinant Rat Cytosolic Fe-S cluster assembly factor NUBP2 (Nubp2)

Product Code	CSB-EP718091RA-B
Abbreviation	Nubp2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q68FS1
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MEAAAEPGNL AGVRHIILVL SGKGGVGKST ISTEALALR HQGKKVGILD VDLCGPSIPH MLHAQ GKAVH QCDSGWVPVF VDQEQSISLM SVGFLENPD EAVVWRGPKK HALIKQFVSD VAWGELDYLVDTPPGTSDE HMATVEALRP YKPLGALVVT TPQAVSIGDV RRELTFCKKT GLQVIGVIEN MSGFACPHCA ECTNVFSSGG GEELARLAGV PFLGSVPLDP QLTRSLEEGR DFIQEFPKST AYSALTSIAH KVLHQMPALC S
Source	E.coli
Target Names	Nubp2
Protein Names	Recommended name: Cytosolic Fe-S cluster assembly factor NUBP2 Alternative name(s): Nucleotide-binding protein 2 Short name= NBP 2
Expression Region	1-271
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.