



Recombinant Rickettsia typhi Elongation factor Ts (tsf)

Product Code	CSB-EP718155RNE
Abbreviation	tsf
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q68XV6
Product Type	Recombinant Protein
Immunogen Species	Rickettsia typhi (strain ATCC VR-144 / Wilmington)
Purity	≥85% (SDS-PAGE)
Sequence	MSEINISAAA VKELREKTGA GMMDCKKALI ETSGNFEEAI DFLRKKGLAA AVKKFGRIAS EGLTAIKING LTSVVIEVNS ETDFVARNEQ FQNLVKDIVN LAIIAQNIDA LKISKMQSGK SVEEEIENI ATIGENLTLR RMDILEISNG AIGSYVHNEV VPHLGKISVL VGLESNAKDK AKLAALAKQI AVHVAGNNPQ SIDTSLDQA LVERERKVFF EKSKEEGKPD HIEKMVEGR IRKFFSEVVL LQQNFLFEPK LTVAEVIKNA EQELGAEIKI TKFIRYALGE GIEHEEKNA DEVASITKC
Source	E.coli
Target Names	tsf
Protein Names	Recommended name: Elongation factor Ts Short name= EF-Ts
Expression Region	1-309
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.