



Recombinant Gorilla gorilla gorilla Probable DNA dC->dU-editing enzyme APOBEC-3C (APOBEC3C)

Product Code	CSB-MP718165GGZ
Abbreviation	APOBEC3C
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q694B5
Product Type	Recombinant Protein
Immunogen Species	Gorilla gorilla gorilla (Western lowland gorilla)
Purity	≥85% (SDS-PAGE)
Sequence	MNPQIRNPMK AMYPGTFYFQ FKNLWEANDR NETWLCFTVE GIKRRSVVSW KTGVFRNQVD SETHCHAERC FLSWFCDDIL SPNTNYQVTW YTSWSPCPEC AGEVAEFLAR HSNVNLTIPT ARLYYFQDTD YQEGRLSLSQ EGVAVKIMDY KDFKYCWENF VYNDDEPFKP WKGLKYNFRF LKRRRLQEILE
Source	Mammalian cell
Target Names	APOBEC3C
Protein Names	Recommended name: Probable DNA dC->dU-editing enzyme APOBEC-3C EC= 3.5.4.-
Expression Region	1-190
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.