



Recombinant *Oryza sativa* subsp. japonica Putative zinc finger CCCH domain-containing protein 58 (Os09g0305900, LOC_Os09g13530)

Product Code	CSB-MP718179OFG
Abbreviation	Os09g0305900, LOC_Os09g13530
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q69KP0
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MDSAAAERPS SYAVNFPPLL PAPAPAVAGA MGVANHKSVL CMKWREGRCH NGVACRYAHG EEDQRIVPEM RVGGGGTSMH ARSSPPRDGA SSGSTASIAM AACRIEEQRH GRGGESFILP RSRKRQALG SGSARSTAPT PPRAHTTPPC SRSVAAPRAS TSPLRPSPVP PPATLHSAAD VQRSVARALE DFEQRESSSS VFPLAIDIVA EDAMTATSEP SATSDDDAIT TTTSSSTTDA DELDAAVAAP PK
Source	Mammalian cell
Target Names	Os09g0305900
Protein Names	Recommended name: Putative zinc finger CCCH domain-containing protein 58 Short name= OsC3H58
Expression Region	1-252
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.