



# Recombinant *Oryza sativa* subsp. *japonica* Endoglucanase 23 (GLU12)

<b>Product Code</b>	CSB-YP718180OFG
<b>Abbreviation</b>	GLU12
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q69NF5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	A PSAAAAAAGH DYGDALAKSI LFFEGQRSGR LPAAGQRAAW RGDSAVSDGG AAGVDLEGGY YDAGDNVKFG FPMAFTATML AWGVVEFGDA MPPAERAHAA DAVRWATDYL LKTISHPGVI FIQVGDPTKD HGCWERPEDM DTARTVYNIS AARPGSDVAG ETAAALAAAS MVFRDDDPAY AARLLAGARS AFEFADEHKG AYSDDPELRA GGCPFYCDFD GYQDELLWGA AWLRRASKEG TYLDYIQNNG KTLGAEDSTN EFGWDNKHAG INVLVSKEFI DGEVLSLQSY KEFADGFICT LIPESSPHI TYTPGGMIYK PGGSNMQHVT SISFLLLYA KYLSNSSRTV NCGNVSVGPA TLQQLARKQA DYILGDNPMK MSYMGYGDY YPQRIHHRGS SLPSIKSHPQ RIACNDGTPY YNSSSPNPNP LIGAVVGGPG EDDVYEDDRA DFRKSEPTY INAPLVGVLA YLVGNPDPGG GHVRH
<b>Source</b>	Yeast
<b>Target Names</b>	GLU12
<b>Protein Names</b>	Recommended name: Endoglucanase 23 EC= 3.2.1.4 Alternative name(s): Endo-1,4-beta glucanase 23 OsGLU12
<b>Expression Region</b>	30-515
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.