



Recombinant Mouse Male-specific lethal 2 homolog (Msl2)

Product Code	CSB-MP718191MO
Abbreviation	Msl2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q69ZF8
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MNPVNATALY ISASRLVLNY DPGDPKAFTE INRLLPYFRQ SLSCCVCGHL LQDPIAPTNS TCQHYVCKTC KGKKMMMMPKS CSWCKDYEQF EENKQLSILV NCYKKLCEYI TQTTLARDII EAVDCSSDIL ALLNDGSLFC EETEKPSDSS FTLCLTHSPL PSTSEPTADP QASLSPMSES TLSIAIGSSV INGLPTYNGL SIDRFGINIP SPEHPNTIDV CNTVDIKTED LSDNLPPVCD TVATDLCSTG IDICSFSEDI KPGDSLILLSV EEVLRSLTV SNTVECCPNL QPNLEATVSN GPFLQLSSQS LSHNVFMSTS PALHGLSCTA ATPKVAKLNR KRSRSESDSE KVQPLPISTI IRGPTLGASA PVTVKRESKI SLQPIATVPN GGTPPKISK VLLSTKSMKK SHEHGSKKSH SKSKPGILKK DKAVKEKMPS HHFMPGSPTK TVYKKPQEKK GCKCGRATQN PSVLTCRGQR CPCYSNRKAC LDCICRGCQN SYMANGEKKL EAFVPEKAL EQTRLTLGIN VTSIAVRNAS TSTSVINVTG SPVTTFLAAS THDDKSLDEA IDMRFC
Source	Mammalian cell
Target Names	Msl2
Protein Names	Recommended name: Male-specific lethal 2 homolog Short name= MSL-2 Alternative name(s): Male-specific lethal 2-like 1 Short name= MSL2-like 1 Male-specific lethal-2 homolog 1 RING finger protein 184
Expression Region	1-577
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.