



Recombinant *Oryza sativa* subsp. japonica Protein LOL3 (LOL3)

Product Code	CSB-EP718311OFG
Abbreviation	LOL3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6ASS2
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MQSQIVCHGC RSVLRYPSSGA PSVCCALCQA ITTVPPPAPV MEMAHLICGG CRTLLMYTRN ADTVRCSCCS TVNLVRPVNN IAHVSCGQCR TTLMYPYGAP SVKCAICHYI TNTGMNTVAP TPSPMPTSSG SSYNAPPSGS SYNAPPPTSA PTSRPQNVTV VVENPMTVDE KGKLVSNVVV GVTTGK
Source	E.coli
Target Names	LOL3
Protein Names	Recommended name: Protein LOL3 Alternative name(s): Protein LSD ONE LIKE 3 Short name= OsLOL3 Putative zinc finger LOL3
Expression Region	1-186
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.