



# Recombinant Rat 2',5'-phosphodiesterase 12 (Pde12)

<b>Product Code</b>	CSB-YP718323RA
<b>Abbreviation</b>	Pde12
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6AXQ5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MWRLPGRAAL RGVRSVVEQR SPAEATTHEA VRAMERAVVR CVPSEPKLSL SFALADGSHK NMQRDQSEPL GRALSRIATN ALKGHAKVAA AKKSRKNRAH SSGGAACAAT GSEPAATCEP VVKLYYREEA VAEDVLNVDA WQDGAVLQIG DVKYKVERNP PAFTELQLPR YIMAGFPVCP KLGLEFGDPA SSVFRWYKEV KPGAAEPGDS GLASSSHSLQ SSAWIETGVD ERVYTPCNAD IGLRLKLHCT PGNGQRFGPS RELESVCPVE AGPGTCTFDH RHLYTKKVTE DSFIRTVSYN ILADTYAQTE FSRTVLYPYC APYALELDYR QNLIQKELTG YNADLICLQE VDRAVFSDSL VPALEAFGLE GVFRIQHEG LATFYRKSKEF RLLSQHDISF QEALKSDPLH KELLEKLALN PLAQEKVLQR SSVLQISVLQ STTDSSKKIC VANTHLYWHP KGGYIRLIQM AAALVHIRHV SCDLYPGIPV IFCGDFNSTP STGMYHFVIN GSPEDHEDW ASNGEEERCG MSLTHCFKLN SACGEPAYTN YVGGFHGCLD YIFIDLNALE VEQVIPLPSH EEVTTHQALP SVSHPSDHIA LVCDLKWK
<b>Source</b>	Yeast
<b>Target Names</b>	Pde12
<b>Protein Names</b>	Recommended name: 2',5'-phosphodiesterase 12 Short name= 2'-PDE Short name= 2-PDE EC= 3.1.4.-
<b>Expression Region</b>	1-608
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.